

# TRISTEL TRIO WIPES SYSTEM

## Tristel

### Active Constituents:

The active ingredient in the Tristel Wipes System is chlorine dioxide.

### Description:

The Tristel Wipes System is a practical and highly effective way to decontaminate and disinfect non-lumened medical devices, including those that cannot be sterilised by heat. The System provides high-level disinfection in less than 2 minutes.

The Tristel Wipes System is used in hospitals, specialised clinics, GP practices, laboratories, dental practices and vet practices.

The Tristel Wipes System comprises of three easy steps:-

#### - Step One – Pre-Clean

The first step in the decontamination process of instruments is the thorough cleaning of the surface to remove soil and organic matter prior to high-level disinfection. The Tristel Pre-Clean Wipe is impregnated with a low-foaming surfactant system combined with triple enzymes, producing ultra-low surface tension for rapid cleaning of any hard surface.

#### - Step Two – Sporicidal Wipes (High level disinfection step)

The Tristel Sporicidal Wipe incorporates Tristel's patented chlorine dioxide chemistry. It can kill all organisms on a pre-cleaned surface, from which soil and organic matter have been removed, with a contact time of only 30 seconds.

#### - Step Three – Rinse Wipe

The Tristel Rinse Wipe is impregnated with deionised water and a low-level antioxidant which will remove and neutralise chemical residues from a surface that has been decontaminated.

### Directions for Use:

#### Step One - Pre-Clean Wipe:

- a. Put on new gloves.
- b. Dispense a wipe from sachet and wipe over instrument to remove bio burden.
- c. Discard wipe to waste.

#### Step Two - Sporicidal Wipe:

- a. Put on new gloves
- b. Dispense a wipe from sachet.
- c. Unfold the wipe and lay out on the palm of your hand.
- d. If the foam bottle is being used for the first time, depress the pump two to four times to prime the pump. The first output from the foam bottle can be left on the wipe, to be followed by two complete pumps. The foam bottle is then primed for subsequent wipes. For all subsequent wipes, pump two measures of Tristel Activator foam onto the wipe.
- e. Scrunch the wipe for 15 seconds to activate the chlorine dioxide.

## Directions for Use:

### **Step Two – Sporidical Wipe (Continued)**

- f. Wipe the surface of the instrument until it has been covered with Tristel. All areas of the surface must come into contact with the wipe at least once.
- g. Once the entire surface has been wiped and covered with chlorine dioxide, wait 30 seconds.

### **Step Three – Rinse Wipe**

- a. Wipe the instrument to remove any excess foam.

## Presentation:

The Tristel Wipes System is available in two sizes, the Tristel Trio 25 cycle box or the Tristel Trio 50 cycle box.

## Caution:

Do not use if the wipe sachet and/or foam bottle have been damaged.

OTC

# Tristel

